Good Weighing Practice™

The Global Weighing Guideline

Martin Huber
Greifensee
5. Mai 2011

METTLER TOLEDO

GWP®
Good Weighing Practice™

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No written procedure describing calibration of equipment.

For example, thermometers, scales, pH and conductivity meters.

The calibration of instruments and gauges is not done at suitable intervals in accordance with an established written program and with provisions for remedial action in the event accuracy and/or precision limits are not met.

Balances - Qualifications did not include measurement of uncertainty, balance drift, and minimum weight determination, etc.

b. The balance used to weigh more than 20-mg did not comply with the USP 0.1 % requirement for balance measurement uncertainty.

FDA Warning Letter: Risikomanagement

3. The calibration of instruments was not always conducted at suitable intervals.
21 CFR 211.160(b)(4)

[] used for sterilization of []
and a [] used for [] during the qualification of []
[] were not qualified or calibrated. Please provide an impact assessment for
the [] processed in these []
Your response indicates you will qualify and calibrate these pieces of equipment and
identify all critical and non-critical GMP equipment in order to have a more
comprehensive qualification and calibration program.

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FDA Priorities

- Risk Management**
- Better Consumer Information
- Postmarketing Safety Counterterrorism
- Effective Regulation/ Strong Workforce

PHARMACEUTICAL CGMPs FOR THE 21ST CENTURY — **A RISK-BASED APPROACH**

FINAL REPORT

September 2004

FOOD QUALITY AND SAFETY SYSTEMS

A training manual on food hygiene and the Hazard Analysis and Critical Control Point (HACCP) system

ICH HARMONISED TRIPARTITE GUIDELINE

QUALITY RISK MANAGEMENT

Q9

METTLER TOLEDO

FAO/WHO Food Standards

CODEX alimentarius

These General Principles lay a firm foundation for ensuring food hygiene and should be used in conjunction with each specific code of hygienic practice, where appropriate, and the guidelines on microbiological criteria. The document follows the food chain from primary production through to final consumption, highlighting the key hygiene controls at each stage. It recommends a HACCP-based approach wherever possible to enhance food safety as described in *Hazard Analysis and Critical Control Point (HACCP) System and Guidelines for its Application (Annex)*.

14

ISPE GAMP 5

GAMP 5

A Risk-Based Approach to Compliant GxP Computerized systems

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METTLER TOLEDO

This section features a central collage of images and text. At the top left is a smaller version of the lifecycle diagram from the previous slide. Below it is the text 'Routinebetrieb'. In the center is a large blue globe with 'GWP' written on it. To the right of the globe is the text 'Evaluation' and a photo of a man in a white lab coat. Below the globe is the text 'Auswahl' (Selection) and a photo of a laboratory scale. At the bottom of the globe is the text 'Installation' and a photo of a 'Starter-Pac IPac Laboratory Weighing' box. To the left of the globe is the text 'Kalibrierung' (Calibration) and a photo of a technician working at a scale. The METTLER TOLEDO logo is visible on the far left.

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Good Weighing Practice™

METTLER TOLEDO

GWP®
Good Weighing Practice™

GWP® Verification

Customer	Customer Id#	
Contact	Mr. Sample	
Department	Q21 Laboratory	
Zip Code/City	112211 Escondido	
Site	Escondido	Date: 20 July 2010

Instrument Information

GWP® Verification No.	104-20100009-130
Balance Type	AT261
Supplier	
Identification No.	068447
Maximum Capacity	200g
Readability	0.1 mg
Maximum Capacity Fine Range	62 g
Readability Fine Range	0.01 mg
Automatic Internal Adjustment	no
Hazardous Zone	
Use Application	

Requirements/Achievements	Value	Unit
Weighting Accuracy	0.1	%
Minimum Weighing Load Required	200	g
Minimum Weight Required	45	mg
Minimum Weight Achieved	15	mg
Safety Factor Required	1.5	
Safety Factor Achieved	3.0	
Minimum Weigh Certificate	1.0000000000	
Standards & Regulations	USP	
Expansion Factor k	3	

Instrument meets accuracy requirements.
Instrument meets safety factor requirements.

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GWP®
Good Weighing Practice™

Risk Assessment

Weighting Accuracy	0.1 %
Impact of wrong Measurements on Business Process	high
Impact of wrong Measurements on Humans, Animals and on the Environment	medium
Wrong Measurements can be discovered very easily	no
Applicable Standards and Regulations	USP/USMP 923

Recommended Performance Verification Activities

Item	Frequency
Calibration (by user)	before a year
Linearity (by user)	quarterly
Responsibility (by user)	monthly
Sensitivity (by user)	daily

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Recalibration Interval of Weights*

Weight 1	yearly
----------	--------

Test Tolerances*

Weighting	Working Limit	Control Limit
Weight 1	66.7 mg	100 mg
Weight 2	66.7 mg	100 mg
Weight 3	66.7 mg	100 mg

*Note: In case of a change of the control limit, the tolerance limit must be adjusted accordingly.

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